Easy Test For Drinking Water

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SIMPLE TEST FOR DRINKING WATER

Everyone knows that water sustains life but not everyone knows that water can also endanger life as ingestion of or exposure to contaminated water can cause several diseases. Therefore, safe drinking water is essential to good health but conventional methods of water quality testing are dependent on access to sophisticated laboratories which are often not accessible to people of distant areas. Keeping these facts in view this simple micro biological testing kit is designed so that a lay man can easily test microbiological quality of his/ her drinking water.

Components of Microbiological Testing Kit

Following components are required for this test

I	Testing bottle	01
ii	Spirit bottle	01
iii	Bleaching solution	01
iv	Cotton role	01
v	Spirit lamp/ candle	01
vi	Marker	01
vii	Lighter	01
viii	Thermometer	01
ix	Brochure	01

Basic Pattern:

This kit is designed to determine the presence-absence (PA) of bacteria which provides an indication for sanitary quality and level of contamination of water.

Please set the whole kit on safe place near water source.

Using Instructions:

Step-1: Sterilize your hands using spirit from the spirit bottle.

Step-2: Open the tap for few minutes.

Step-3: Spray spirit on tap and flame the tap using lighter.

Step-4: When cool, open the tap again and keep the water flow normal to avoid splashing.

Step-5: Turn "ON" the flame of the spirit lamp.

Step-6: Bring testing bottle out from its box near the flame immediately. The fill the bottle up to the mark and close it tightly quickly to avoid contamination. Shake the bottle well.

Incubation:

Place this bottle in the incubator at 35°C for 24-48 hours or if electricity is unavailable, check the temperature of the room or kitchen using thermometer and place the bottle at suitable place where the temperature is above 25°C for 48 hours.

Observations:

After 24-48 hours check for the color change of the water in the testing bottle and compare to evaluation chart for level of bacterial contamination.

2-1. Pale yellow: no bacterial contamination.

2-2. Black: bacterial contamination.

3. Bottle can be 1/4th, 2/4th, 3/4th or full black depending on the level of bacterial contamination in the water.

Water from Distant Places:

If the water source is a river, stream, pond, well or other places away from home, use the sampling bottles the following way:

- After following the step-1, immediately open the bottle near the flame and fill the sampling bottle and close the tap quickly.
- Before testing, shake the sampling bottle 25 times; sterilize outer surface of bottle use the cotton swab soaked in spirit.
- Near flame pour the water from sampling

to bottle and fill up to mark. Then close the cap quickly and incubate the testing bottles in the incubator (35°C) or at room temperature as mentioned previously.

Disposal of Media:

After using the bottle for testing, use the bleaching solution for disposal of media. Add a few drops of bleaching solution in the used testing bottle, close the cap and mix by shaking and let it stand for 10-15 minutes. Then discard the solution and wash the bottle. At the end wash your hands with soap properly.

Remedial Measures:

If your water is contaminated, boil the water or filter the water before drinking. For the suggestions on how to stop the contamination consult Pakistan Council of Research in Water Resources (PCRWR) at the given telephone number or address.

EVALUATION CHART FOR BACTERIAL CONTAMINATION

OBSERVATION	RESULTS	LEVEL	REMARKS
	→ No Change in Color	→ No bacterial contamination	→ Fit for drinking purpose
	→ Black coloration in 1/4th of bottle	→ Slight contamination	→ Unfit for drinking purpose
	→ Black coloration in 1/2th of bottle	→ Slight higher contamination	→ Unfit for drinking purpose
	→ Black coloration in 3/4th of bottle	→ Fair higher contamination	→ Unfit for drinking purpose
	→ Full black	→ Highly contamination	→ Unfit for drinking purpose